

How to isolate DNA; Large Scale.

Jane Yeadon

Background

It is often necessary to isolate genomic DNA from *Neurospora* strains, to check for particular single nucleotide polymorphisms or to screen transformants for inserted sequences. Both of these can be done using PCR and so require only a small amount of relatively pure DNA. However, if several Southern blots must be done, or a genomic library made of a strain that differs from that sequenced by the Genome project, it may be desirable to isolate a large quantity of good quality DNA

Procedure

Neurospora genomic DNA is prepared using an adaptation of the method of Schectman (1986). A conidial suspension of the desired strain is inoculated into 50 ml of liquid Vogel's medium (supplemented as necessary) in a 250 ml flask and shaken at 25°C for 2-3 days, until the culture is saturated. If the culture is Cot⁺, growth can be at 34°C for a shorter time. We usually keep the cultures in the dark, to discourage conidiation, and give them a good swirl by hand once to twice a day to prevent mycelium growing up the sides of the flask.

The mycelium is harvested by pouring the culture through a funnel lined with two layers of sterile muslin, or passing through a filter paper in a Buchner funnel. Excess liquid is removed by pressing mycelium between layers of paper towel. At this point the mycelial mat can be stored at -70°C. Where possible, the mycelium should be freeze-dried for 24-36 hours, after which it can be stored at -20°C until needed. If using mycelium that had not been freeze-dried, of course you need to use more. The procedure is equally successful, although more difficult, as wet mycelium must be ground under liquid nitrogen.

1g of dry mycelium is ground to a fine powder in a pestle and mortar. The powder is added slowly (to minimise clumping) to 10 ml of *Neurospora* SDS buffer (0.15M NaCl, 0.1M EDTA, 2% SDS, pH 9.5) in a 125 ml flask. Add 0.5 ml 2mg/ml proteinase K solution and incubate at 37°C for up to 24 hrs with gentle agitation.

10 ml sterile water is added to the slurry and cellular debris is removed by centrifugation at 10,000 rpm/10min/4°C. The supernatant is extracted 3-5 times with phenol saturated with Tris-HCl (pH 8.0) and once with water-saturated chloroform.

The DNA is precipitated by the addition of 2.5 volumes of absolute ethanol and removed with a sterile glass hook to an Eppendorf tube. If this doesn't work for you, you can centrifuge briefly and remove the supernatant. After air-drying briefly, add 1ml sterile 10mM Tris-EDTA, pH 8.0 and 30µg. Do not use any vigorous means of resuspending the DNA, as it will shear easily. After incubation at 37°C for several hours (or overnight) with occasional mixing by inversion, the DNA should dissolve by itself. To remove the RNase, add 80µg of proteinase K, and continue incubation at 37°C for a further 2 hours.

The opalescent liquid is then extracted 3 times with phenol/Tris and twice with chloroform/H₂O. Precipitate DNA with 2.5 volumes of ethanol (brief centrifugation only required), wash with 70% ethanol, dry briefly and resuspend in up to 2 ml 1xTE.

This method has yielded 1-2 ml of genomic DNA sufficiently concentrated that 2µl are sufficient for a non-radioactive blot. DNA I have made as long ago as 1992 is still in very good condition after extended storage at 4°C.